
Implants + Instruments

# Distraction Devices

Enorale Distraktoren für die Ober- und Unterkiefer  
Intraoral Distraction Devices for Upper and Lower Jaw  
Distractores intraorales para la mandíbula inferior y superior  
Distracteurs intra-oraux pour mâchoire inférieure et supérieure  
Distrattori intraorali per il mascellare inferiore e superiore

**medicon**<sup>®</sup>

Not all products illustrated in this documentation are available for sale in all countries.  
Not for all declared products in this documentation Medicon is the legal manufacturer.  
Please contact your local distributor for more information concerning availability and regulatory information (e.g. classification) of the mentioned products.

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
Distraction devices for intraoral application generally facilitate the lengthening of an under-developed mandible on one or both sides in the horizontal and ascending mandibular ramus. Special distraction devices can be used to obtain bone even from the median section of the mandible. In each case, transantral distractors are arranged in pairs parallel to each other. The particular design of the miniplates for osteosynthesis allows the system to be highly versatile for its application to different bony surfaces; as well as to achieve an optimum "distraction vector" for implant prosthetic rehabilitation purpose.

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

The horizontal distraction device can be used to lengthen the horizontal mandibular ramus. The horizontal distraction device comprises a distraction cylinder with a T-shaped mini-osteosynthesis plate for monocortical attachment of the distraction device medially to the mandibular transverse osteotomy and a distraction pole (that can be lengthened by turning the screws) with an L-shaped mini-osteosynthesis plate for mono- or bicortical attachment of the distraction device distally to the osteotomy. The distraction pole is pushed out of the cylinder by turning the screw at the end of the distraction cylinder to the left. This enables the mandible to be lengthened as desired. If the horizontal mandibular ramus is to be lengthened, the horizontal distraction device should be fixed parallel to the occlusion plane with monocortical screws and the transbuccal drill and screw set (68.76.90) or angular screwdriver (68.30.01). Structures such as tooth germs, tooth roots or nerves can be preserved through mono-cortical attachment of the device. If the ascending ramus is to be lengthened, then the distraction device should be positioned obliquely to the occlusion plane such that non-occlusion occurs as a result of mandibular lengthening. Particular attention should be paid to the position of the distraction cylinder in this instance so that the distraction device can still be used with the screwdriver.

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---

**Indications:** Hemifacial microsomia possibly with airway obstruction  
Hemifacial microsomia possibly with facial asymmetries  
Pierre-Robin-Syndrome • Treacher-Collins-Syndrome • Crouzon-Syndrome  
Extreme facial asymmetries, shortened collum mandibulae due to avascular necrosis or post-traumatic  
Mandibular retrognathia possibly with impairment of the mandibular joint, arthrosis deformans of the mandibular joint

**Goal:** Horizontal progenia of the mandibular ramus for aperture of the airways and for primary mandibular shifts

**Later measures:** Vertical distraction for the creation of a mandibular angle

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
Intraoral distraction up to 20 mm in the milk or changing dentition for lengthening of the horizontal mandibular ramus. Application of the distraction device parallel to the occlusion plane [fig. 1]. Two months after completion of the distraction of 8 mm with ossification of the distraction distance [fig. 2]. Intraoral distraction up to 25 mm in the final dentition to lengthen the horizontal dysplastic mandibular ramus. Application of the distraction device parallel to the occlusion plane [fig. 3]. Two months after completion of the distraction of 25 mm with ossification of the distraction distance [fig. 4]. Intraoral distraction up to 20 mm or 25 mm in the milk, changing or final dentition to lengthen the ascending mandibular ramus.

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


Application of the distractor parallel to the direction of the distraction or to the growth direction of the dysplastic mandibular ramus [fig. 5]. A mandibular progenia with non-occlusion is achieved [fig. 6]. Pre-bending of the L- and T-shaped mini-plates on the horizontal distraction device while simultaneously opening of the distraction cylinder by 5 mm in order to control the transversal dimension of the deformed mandibular ramus [fig. 7]. Extension of the osteotomated, deformed and shortened mandibular ramus by means of pre-bending of the mini-plates of the distraction device with formation of a wedge-shaped osteotomy cleft [fig. 8]. Subsequent distraction of the deformed mandibular ramus in horizontal direction to adjust the transversal dimension and the mandibular middle line [fig. 9].

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


The vertical distraction device capable of lengthening the ascending ramus in various directions also consists of a distraction cylinder to which an L-shaped mini osteosynthesis plate is fixed caudally to the osteotomy for monocortical attachment of the distraction device, and a concealed distraction pole (that can be lengthened by turning the screws) with a T-shaped mini-osteosynthesis plate for mono- or bicortical attachment cranially to the osteotomy. If the ascending mandibular ramus is to be lengthened, then the vertical distraction device must be accurately positioned in the direction of the distraction. The screws can be fixed bicortically above the Nervus alveolaris inferior, but preferably monocortically in the vicinity of the nerve. A variety of vertical distraction directions can be adjusted with this device so that progenia of the mandibular corpus is no longer essential. Only the effect of non-occlusion is apparent. Unrestricted use of these vertical distraction devices is feasible by changing the drive direction with a cardan gear and universal joint. Care should be taken to ensure that the isolated operating rod in the vestibulum does not incline by more than 30° to the cardan gear on turning to the left. A distraction device can be used in children from 4 years of age. The vertical 15 mm distraction device is particularly suitable in this case. The vertical 20 mm distraction device is used in older patients.

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---

**Indications:** Hemifacial microsomia  
Shortened ascending mandibular ramus  
Open bite based on the shortened ascending mandibular ramus  
Pierre-Robin-Syndrome • Treacher-Collins-Syndrome  
Craniofacial microsomia

**Goal:** Vertical distraction to lengthen the ascending mandibular ramus, creation of a prominent mandibular angle

**Later measures:** Orthodontic after-treatment

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Intraoral distraction up to 20 mm in the milk- and changing dentition to lengthen the ascending mandibular ramus. The position of the distraction cylinder plots the distraction direction. Two months after completion of the distraction, ossification of the distraction distance can be seen. The vertical 20 mm distraction device is used on the ascending ramus in the event of minor deformity and when there is sufficient space for distractor adaptation. A bone adaption surface of 10×25 mm is required in order to attach the distraction device. The vertical 15 mm distraction device is used on the ascending ramus in the event of extensive deformity and in confined spaces for distractor adaption. A bone adaption area of 10 ×25 mm is required in order to attach the distraction device [fig. 1-3].

▶ Simple activation of the distraction device by means of a screwdriver: 1 counter-clockwise rotation = 0.4 mm [fig. 4].

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




When performing a distraction osteogenesis an exact planning of each step of the surgery is essential. Part of this planning is the calculation and application of the angle of the distraction osteogenesis. The combination of the positioning rod with a cotter pin offers the advantage of a very easy adoption of the required distraction angle during the surgery.

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

The class II distraction device is used to lengthen an excessively short mandible in adults. This distraction device with its attached, reusable positioning device is constructed so as to facilitate parallel mandibular progenesis without altering the position of the condyle mandibular. The cephalo position of the distraction device, planned on lateral X-ray, can be transferred during surgery with the aid of the positioning device and its individual, occlusal, synthetic splint.

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**Advantages**

Reduced risk of alveolar nerve damage  
Gradual mandibular advancement with slow stretching of soft tissue  
Dynamic occlusal adjustment  
Occlusion control and adjustment feasible by orthodontics  
Prevention of an immediate surgical relapse  
Reduced risk of TMJ unit problems

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
---

**Indications:** Mandibular retrognathia • Angle Class II deep bite

**Contra-indications:** Open bite • Facial asymmetry • Laterognathia

**Prerequisites:** Orthodontic pre-treatment • Congruent dental arch  
Deep bite or parallelism between maxilla and mandibular occlusion plane • Middle line displacement < 3 mm

**Goal:** The aim of this distraction technique is the gradual advancement of the mandible in an eugnathic position with functional loading of the mandible

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
Horizontal lengthening of the mandibular corpus is feasible with this distraction device. It is preferably used in adults for the correction of mandibular retrognathies following orthodontic pre-treatment [fig. 1]. This device is similar to the horizontal distraction device but has an angular, distal, L-shaped mini-osteosynthesis plate for the distal attachment of the device to the exterior of the mandibular ramus. This angular mini-osteosynthesis plate facilitates the anatomically correct positioning of the device in the region of the anterior margin of the ramus so that the distraction cylinder is fixed flush to the dental arch, approximately 1 cm below the occlusion plane [fig. 2-3]. The distraction cylinder is thus placed in a parallel position, preventing lateral displacement of the head of the jaw [fig. 4].

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

A positioning device is used for a more accurate parallel positioning of the distraction cylinder (69.00.65). The front ends of the left and right distraction cylinder are fixed in this positioning device by means of imbus screws. On the lateral distal X-ray, the horizontal position of the distraction device is parallel to the occlusion plane. This position is transferred to the mandibular model using an individual, occlusal, synthetic splint. The splint is temporarily fixed to the bracket with wire during surgery. The level plane of the splint corresponds to the distraction direction. The positioning device with its occlusal plate is screwed onto the splint so that the position of the distraction device confirmed by distal X-ray and in the model is transferred exactly above the splint and the positioning device at the operating site [fig. 5-6]. After removal of the positioner distraction is effected under orthodontic control until a correct bite is reached [fig. 7-8].

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





The median mandibular distraction device can be used for the transversal widening of the dental arch up to 10 mm. It comprises a closed distraction cylinder with L-shaped mini-osteosynthesis plates for bicortical attachment of the distraction device laterally to the vertical, interdental or median mandibular osteotomy, and a distraction pole (that can be lengthened by turning the screw) also with L-shaped mini-osteosynthesis plates for bicortical attachment of the distraction device to the other side of the vertical osteotomy. Median, transversal mandibular distraction is indicated for crowded front teeth. Orthodontic practice to date has advocated premolar extraction. The ensuing gaps were preferably closed by resolving the crowding but, in many cases, also by additional retrusion of the front of the mandible, once again reducing upper lip prominence. Such tooth extractions are now no longer necessary. The crowded teeth in the median gaps caused by distraction slowly move closer together. The distraction distance should be equivalent to the amount of space calculated orthodontically to resolve the crowding. The advantage of this treatment is that protrusion of the front dentition is preserved, which is essential for upper lip prominence. In future, premolar extractions on orthodontic grounds should only be rarely indicated.

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
Median mandibular distraction is indicated for transversal narrowing of the dental arch and crowding of the front teeth [fig. 1]. Bicortical attachment of the distraction device with mini plates on stability reasons [fig. 2]. Daily distraction distance 2-3×0.4 mm equivalent to 2-3×360° counter-clockwise rotations of the screwdriver as from the 3rd or 4th day after surgery. Slow expansion of the osteotomy fissure [fig. 3], which is completely ossified after 3 months. The mandibular front teeth are subsequently moved orthodontically into the gaps to resolve the crowding and to form the dental arch [fig. 4]. The maxillary dental arch is also to be formed congruently to the mandible.

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
The distraction device is operated by turning the screw to the left. A complete rotation corresponds to distraction of 0.4 mm. Daily distraction distances of 3×0.4 mm in children and 2×0.4 mm in adults should be achieved.

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
Transantral distraction of the maxilla allows gradual advancement of the maxilla by up to 25 mm without letting the distractors become visible. In each case, transantral distractors are arranged in pairs parallel to each other. They are activated with screwdrivers using two jointed extension pieces in the frontal vestibule of the maxilla. These extension pieces are removed at the end of the distraction so that the distractors are completely covered by mucosa until they are removed.

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**Indications:** Patients with congenital and acquired maxillary retrognathia, for example, in cheilognathouranoschisis, hypoplasia of the maxilla in syndromes or following severe midfacial comminuted fractures.

**Condition:** Bony continuity of the piriform aperture, the canine fossa and the posterior wall of the maxillary sinus is required as the distractors are fixed to these points. Moreover, the Le Fort I osteotomy must be modified and the osteotomy performed in front of the maxillary tuberosity so that the bony pterygoid process is not injured. Stability of the posterior bony wall of the sinus is a further requirement for correct fixation of the transantral distractor. If these conditions are not met, for example in the case of cheilognathouranoschisis patients, reconstruction of these sites by means of bone grafts must be performed prior to the distraction.

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
The tips of the distractor pins are first pressed into the bony posterior sinus wall lateral to the neurovascular bundles. The distractors are then extended with the screwdriver and adjusted to the sagittal length of the maxilla. The straight miniplates at the front edge of the distraction cylinder are each fixed to the piriform aperture with three miniscrews. Both distractors are positioned stably parallel to each other. Depending on the distraction direction previously determined on the lateral radiograph, the angled miniplates of the two distractors are fixed to the canine fossa of the maxilla in such a way that only a narrow osteotomy gap remains at the Le Fort I level. The maxilla can thus be advanced parallel to the osteotomy line. However, it is also possible to fix the mobile maxilla further caudally so that a large gap arises at the Le Fort I level. In this way, distraction can be performed in the anterior caudal direction.

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The jointed extension piece allows problemfree distraction of the maxilla through the maxillary vestibule at an angle of 70°. At the end of the distraction when the maxilla has been advanced, normal occlusion has been obtained and after intraoral removal of the extension piece, the distractor is completely covered by soft tissue and can remain in place (crossing the lumen of the sinus without irritation) for an unlimited period. This rules out recurrence due to scarring. The two removable extension pieces, which are fixed intraoperatively to the front end of the distraction cylinder, penetrate the vestibular mucosa and so make it possible to activate the distractor intraorally. They are secured in the maxillary vestibule with a ring and a nonabsorbable suture, which is fixed in the mucosa. At the end of the distraction process, the two extension pieces are removed. Thus the distractor is completely surrounded by tissue and is no longer in communication with the oral cavity.

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
Applied intraorally, this distraction device is used to reconstruct one or more sectors of the mandible with a vertical bone deficit, for implant-prosthetic purposes. Made of titanium, the distractor consists of a distraction rod with microplates (0.7 mm thick) for monocortical fixation to the basal mandible, a distraction cylinder with microplates (0.7 mm thick) for monocortical fixation to the mobile bone segment and a screwdriver to activate the distraction mechanism. The distractor is lengthened and the mobile bone is consequently elevated when the screw is turned: one full turn of the screw (360°) in the anti-clockwise direction lengthens the distractor by 0.3 mm, maximum lengthening distance being 13 mm. Distraction rate: 0.9 mm per day (may be reduced in case of clinical evidence of impairment of the distracted soft tissues). Distraction steps: 1 or 2 steps per day. The telescopic construction with a fully enclosed distraction mechanism reduces to a minimum the interference with the callus formation inside the distraction gap and with the healing of the surrounding soft tissue.

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Horizontale Dimension: ~15,0 mm • Vertikale Dimension: ~8,0 mm • Transversale Dimension: ~4,0-5,0 mm

---

**Indications:** Dento-alveolar mandibular deficit:

Sectorial mandibular resection for benign or malignant neof ormation  
As a result of trauma

Mandibular atrophy:

Frontal sector • Rear sector • Total

The basal bone segment where to apply the device must fulfill at least following demands:

Horizontal dimension: ~15.0 mm • Vertical dimension: ~8.0 mm • Transversal dimension: ~4.0-5.0 mm

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Dimension horizontale: ~15,0 mm • Dimension verticale: ~8,0 mm • Dimension transversale: ~4,0-5,0 mm

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Application of the distractor. Post-operative latency period to distraction 5-7 days, depending on the age of the patient [fig. 1]. Distraction rate is 0.9 mm per day (depending on the clinical condition of the soft tissues subjected to expansion). Distraction frequency: 1-2 times a day [fig. 2]. Stabilization period: 2-3 months. Stabilisation period: the correct position of the distractor corresponds with the mandibular arch [fig. 3]. After removal of the distractor the expanded part will be completely ossified. The dental implants can be inserted immediately after removal of the distractor in the same procedure [fig. 4].

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



The channel enables a secure guidance of drill and screwdriver and supports screwing in areas with difficult access.

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

In all situations in which drilling and screwing operations in the mandibula are hampered by the cheek, this instrument set is fast and simple to use. It is small and convenient, as all impeding handles can be unscrewed if required. It can be left in place and used repeatedly as surgery progresses, its compact size prevents it from impeding the work of the surgeon. A further advantage is that the cheek remains freely movable in all directions although the plug-in sleeve is applied. This allows an extensive drilling and screwing in the lateral area of the mandibula.

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Angular screwdriver for plate osteosynthesis in not easily accessible regions of the mandible and mid-face instead of a transbuccal/transfacial softtissue perforation directly through the intraoral supra-/infraorbital or coronal approach to a fracture.

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—MEDICON—  
**implantArt**  
 —OSTEOSYNTHESIS—

**MICRO**  
**1.5**

a	"EASY GRIP"		"MASTER GRIP"	
	Self-tapping Ø 1,5 mm	Emergency screw Ø 1,9 mm	Self-tapping Ø 1,5 mm	Emergency screw Ø 1,9 mm
3,5 mm	<b>68.83.53A</b>		<b>68.83.03A</b>	
5,0 mm	<b>68.83.55A</b>		<b>68.83.05A</b>	
6,0 mm	<b>68.83.56A</b>	<b>68.83.86A</b>	<b>68.83.06A</b>	<b>68.83.36A</b>
7,0 mm	<b>68.83.57A</b>		<b>68.83.07A</b>	

a	"POWER GRIP"		"STAR GRIP"	
	Self-tapping Ø 1,5 mm	Emergency screw Ø 1,9 mm	Self-tapping Ø 1,5 mm	Emergency screw Ø 1,9 mm
3,5 mm	<b>68.85.43A</b>		<b>68.85.03A</b>	
5,0 mm	<b>68.85.45A</b>		<b>68.85.05A</b>	
6,0 mm	<b>68.85.46A</b>	<b>68.85.66A</b>	<b>68.85.06A</b>	<b>68.85.26A</b>
7,0 mm	<b>68.85.47A</b>		<b>68.85.07A</b>	


Verpackungseinheit 10 Schrauben im Magazin (nicht steril)  
 Packing unit 10 screws in storage tray (non-sterile)  
 Paquete de 10 tornillos en cartucho (no estéril)  
 Emballage unitaire de 10 vis en cartouche (non stériles)  
 Confezione da 10 viti nella cartuccia (non sterile)

STRYKER		Ø 1,2 mm
	a = 6 mm    b = 48 mm	<b>68.76.32</b>
	a = 8 mm    b = 50 mm	<b>68.76.33</b>
	a = 18 mm    b = 50 mm	<b>68.76.34</b>
DENTAL		Ø 1,2 mm
	a = 18 mm    b = 50 mm	<b>68.76.35</b>

Selbstschneidend  
**Self-tapping**  
 Autoroscante  
 Autotaraudeuse  
 Autofilettante


Notfallschraube (selbstschneidend)  
**Emergency screw (self-tapping)**  
 Tornillo de emergencia (autoroscante)  
 Vis d'urgence (autotaraudeuse)  
 Vite d'emergenza (autofilettante)


**Schraubendreher und Einsätze für Winkelschraubendreher**  
**Screwdrivers and Inserts for Angular Screwdriver**  
**Destornilladores e insertos para destornilladores de ángulo**  
**Tournevis et embouts pour le tournevis angulaire**  
**Cacciaviti ed inserti per cacciaviti angolati**

**MICRO**  
**1.5**


Schraubendreherklingen  
 Screwdriver blades  
 Hojas de atornillador  
 Lames de tournevis  
 Lame per cacciavite


Schraubendrehergriffe  
 Screwdriver handles  
 Mangos para atornillador  
 Manches de tournevis  
 Manici per cacciavite

Schraub-Einsätze für Winkelschraubendreher  
 Screwing inserts for angular screwdrivers  
 Insertos para destornillador para destornilladores de ángulo  
 Embouts pour le tournevis angulaire  
 Inserti per strumenti per viti per cacciaviti angolati

Vorbohrer-Einsatz für Winkelschraubendreher  
 Drill insert for angular screwdrivers  
 Inserto de broca para destornilladores de ángulo  
 Embout perforateur pour tournevis angulaire  
 Perforatore per cacciaviti angolati


**Titan-Schrauben und Spiralbohrer**  
**Titanium Screws and Twist Drills**  
**Tornillos de titanio y brocas espirales**  
**Vis en titane et forets hélicoïdaux**  
**Viti in titanio e perforatori elicoidali**



—MEDICON—  
**implantArt**  
 —OSTEOSYNTHESIS—

**MINI**  
**2.0**



<b>a</b>
5,5 mm
7,0 mm
9,0 mm

“EASY GRIP”	
	
Self-tapping	Emergency screw
Ø 2,0 mm	Ø 2,3 mm
<b>68.75.25A</b>	<b>68.75.42A</b>
<b>68.75.27A</b>	<b>68.75.43A</b>
<b>68.75.29A</b>	

“MASTER GRIP”	
	
Self-tapping	Emergency screw
Ø 2,0 mm	Ø 2,3 mm
<b>68.75.05A</b>	<b>68.75.22A</b>
<b>68.75.07A</b>	<b>68.75.23A</b>
<b>68.75.09A</b>	

<b>a</b>
5,5 mm
7,0 mm
9,0 mm


“POWER GRIP”	
	
Self-tapping	Emergency screw
Ø 2,0 mm	Ø 2,3 mm
<b>68.94.01A</b>	<b>68.94.09A</b>
<b>68.94.02A</b>	<b>68.94.10A</b>
<b>68.94.03A</b>	


“STAR GRIP”	
	
Self-tapping	Emergency screw
Ø 2,0 mm	Ø 2,3 mm
<b>68.78.05A</b>	<b>68.78.22A</b>
<b>68.78.07A</b>	<b>68.78.23A</b>
<b>68.78.09A</b>	

Verpackungseinheit 10 Schrauben im Magazin (nicht steril)  
 Packing unit 10 screws in storage tray (non-sterile)  
 Paquete de 10 tornillos en cartucho (no estéril)  
 Emballage unitaire de 10 vis en cartouche (non stériles)  
 Confezione da 10 viti nella cartuccia (non sterile)


		
STRYKER		Ø 1,45 mm
	a = 11 mm    b = 44 mm	<b>68.76.73</b>
	a = 18 mm    b = 70 mm	<b>68.76.63</b>
	a = 18 mm    b = 105 mm	<b>68.76.20</b>

		
DENTAL		Ø 1,45 mm
	a = 11 mm    b = 44 mm	<b>68.79.89</b>

Selbstschneidend  
**Self-tapping**  
 Autoroscante  
 Autotaraudeuse  
 Autofilettante


Notfallschraube (selbstschneidend)  
**Emergency screw (self-tapping)**  
 Tornillo de emergencia (autoroscante)  
 Vis d'urgence (autotaraudeuse)  
 Vite d'emergenza (autofilettante)





**Schraubendreher und Einsätze für Winkelschraubendreher**  
**Screwdrivers and Inserts for Angular Screwdriver**  
**Destornilladores e insertos para destornilladores de ángulo**  
**Tournevis et embouts pour le tournevis angulaire**  
**Cacciaviti ed inserti per cacciaviti angolari**

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© Copyright 12/2024, MEDICON eG, Tuttlingen, Germany